

## UNIT 3

### Introduction to Precipitation Titrations

Precipitation titrations involve the **formation of an insoluble precipitate** during a chemical reaction between the analyte and the titrant. These are mainly used for the **determination of halides ( $\text{Cl}^-$ ,  $\text{Br}^-$ ) and certain metal ions** using **silver nitrate ( $\text{AgNO}_3$ )** as a titrant. The common indicators used in such titrations depend on the nature of the precipitate.


### Mohr's Method

#### Definition:


Mohr's method is a **direct precipitation titration** used for the **determination of chloride and bromide ions** by titration with **standard silver nitrate solution**, using **potassium chromate ( $\text{K}_2\text{CrO}_4$ )** as an indicator.

#### Principle

- When **silver nitrate** is added to a solution containing **chloride ions**, **white precipitate of silver chloride ( $\text{AgCl}$ )** forms:



- After all chloride ions are precipitated, **excess silver ions** react with **chromate indicator** to form a **red precipitate of silver chromate ( $\text{Ag}_2\text{CrO}_4$ )**, indicating the end point:



#### Reagents

- Titrant:** Standard **silver nitrate ( $\text{AgNO}_3$ )** solution
- Indicator:** **Potassium chromate ( $\text{K}_2\text{CrO}_4$ )** (1% solution)
- Analyte:** Chloride-containing sample (e.g.,  $\text{NaCl}$ , tap water)

#### Procedure

- Pipette a measured volume of chloride solution into a conical flask.
- Add 1–2 mL of **1% potassium chromate indicator**.
- Titrate with **0.1 N  $\text{AgNO}_3$**  from a burette.
- During titration, a **white precipitate** of  $\text{AgCl}$  forms.
- The **end point** is indicated by the **first permanent reddish-brown color** due to the formation of  **$\text{Ag}_2\text{CrO}_4$** .

### Conditions for Accurate Results

- The solution should be **neutral to slightly alkaline** (pH ~7–9); in acidic media, chromate converts to dichromate and in basic media silver hydroxide may precipitate.
- Avoid strong light and temperature changes as silver salts are photosensitive.

### Calculations

From the volume of  $\text{AgNO}_3$  used and its normality, the amount of  $\text{Cl}^-$  present can be calculated using:

$$1 \text{ mL of } 0.1 \text{ N AgNO}_3 \equiv 0.003545 \text{ g}$$

### Applications

- Estimation of **chloride in drinking water, urine, or pharmaceutical preparations**
- Quality control of **sodium chloride injections**
- Analysis of **brine solutions**


## 3. Volhard's Method

### Definition:


Volhard's method is an **indirect precipitation titration** used to estimate **halide ions (e.g.,  $\text{Cl}^-$ ,  $\text{Br}^-$ ,  $\text{I}^-$ )** by **back titration** with **standard thiocyanate solution** in **acidic medium**, using **ferric ammonium sulfate** as an indicator.

### Principle:


- A known excess of **standard silver nitrate ( $\text{AgNO}_3$ )** is added to the halide-containing solution to precipitate all halide ions as  **$\text{AgCl}$ ,  $\text{AgBr}$ , or  $\text{AgI}$** :



- The excess unreacted  **$\text{AgNO}_3$**  is back-titrated with **standard ammonium thiocyanate ( $\text{NH}_4\text{SCN}$ )**:



- At the endpoint, **ferric ions ( $\text{Fe}^{3+}$ )** react with the first excess of thiocyanate to form a **blood-red complex**:



### Medium:

- Acidic (usually nitric acid) to prevent precipitation of  $\text{Fe}(\text{OH})_3$  and to dissolve any formed  $\text{Ag}_2\text{CrO}_4$

**Indicator:**

- Ferric ammonium sulfate (gives red color with  $\text{SCN}^-$ )

**Applications:**


- Estimation of **chloride, bromide, iodide** in **insoluble salts**
- Also used for **silver** determination (by direct titration with  $\text{SCN}^-$ )

**4' Modified Volhard's Method****Definition:**


Modified Volhard's method is an **indirect titration method** used for the estimation of **halide ions ( $\text{Cl}^-$ ,  $\text{Br}^-$ ,  $\text{I}^-$ )** where the precipitated **silver halide** ( $\text{AgCl}$ ,  $\text{AgBr}$ ,  $\text{AgI}$ ) is **filtered off** before the back titration. This modification is made to avoid **interference** caused by the precipitate during endpoint detection.

**Principle:**


1. A known excess of **standard silver nitrate ( $\text{AgNO}_3$ )** is added to the halide-containing solution to form a **precipitate of silver halide**:



2. The precipitate ( $\text{AgCl}$ ) is **filtered off**, and the **filtrate** (which contains the excess unreacted  $\text{Ag}^+$ ) is collected.
3. This excess silver nitrate is then titrated with **standard ammonium thiocyanate ( $\text{NH}_4\text{SCN}$ )** solution:



4. The **endpoint** is detected by the formation of a **blood-red complex** between the first excess of thiocyanate and added **ferric ions**:



---

**Reagents Required:**

- **Standard  $\text{AgNO}_3$  solution** (e.g., 0.1 N)
- **Standard  $\text{NH}_4\text{SCN}$  solution** (e.g., 0.1 N)
- **Ferric ammonium sulfate** ( $\text{Fe}^{3+}$  indicator)
- **Nitric acid** (to maintain acidic medium)
- **Halide-containing sample** (e.g.,  $\text{NaCl}$  solution)

### Procedure:

1. Take the halide-containing sample and add a **known excess of standard AgNO<sub>3</sub> solution** in acidic medium (HNO<sub>3</sub>).
2. Allow complete precipitation of AgCl.
3. **Filter the precipitate** to remove AgCl and collect the **filtrate**.
4. Add **1–2 mL of ferric ammonium sulfate indicator** to the filtrate.
5. Titrate with **standard NH<sub>4</sub>SCN** until a **permanent blood-red color** appears, indicating the endpoint.

### Advantages of Modified Volhard's Method:

- **More accurate** than Mohr's method when interfering substances or colored solutions are present.
- **Filtration step** removes surface-adsorbed SCN<sup>-</sup> or Fe<sup>3+</sup> from AgCl, eliminating endpoint errors.
- Suitable for **colored, turbid, or complex matrices** where direct titration is difficult.

### Applications:

- Estimation of **chloride** in pharmaceutical and biological samples.
- Suitable for samples where **silver halide precipitate** may adsorb indicator or SCN<sup>-</sup>, causing inaccurate endpoints.

### Limitations:

- **Filtration step** adds complexity and time.
- Requires **strict control** of pH and cleanliness to avoid silver ion loss.


## Fajans Method

### Definition:

Fajans method is a **precipitation titration technique** that uses **adsorption indicators** to detect the **end point** of halide titration with **silver nitrate (AgNO<sub>3</sub>)**. It is used primarily for the **direct estimation of halides** like chloride (Cl<sup>-</sup>), using indicators such as **eosin, fluorescein, or dichlorofluorescein**.

### Principle:

- When AgNO<sub>3</sub> is added to a solution containing Cl<sup>-</sup>, a **white precipitate of silver chloride (AgCl)** forms:



- Near the end point, after all Cl<sup>-</sup> is precipitated, **excess Ag<sup>+</sup>** is adsorbed on the surface of the AgCl precipitate, making it **positively charged**.

- The **negatively charged indicator anion** (e.g., fluorescein<sup>-</sup>) is then **adsorbed onto the AgCl surface**, forming a **pink or red layer**, signaling the **end point**.

#### Key Conditions:

- **pH** should be around **neutral (6–8)** for proper functioning of the indicator.
- **Indicator** must be **adsorbed only after the equivalence point**.
- The titration should be **performed in dim light** to prevent photodecomposition of AgCl.

#### Indicators Used:

- **Fluorescein, eosin, or dichlorofluorescein**
- These give a visible **color change** upon adsorption onto AgCl

## Estimation of Sodium Chloride


#### Reagents Required:

- Standard **0.1 N AgNO<sub>3</sub>**
- Sample solution containing **NaCl**
- **Eosin Y** or **fluorescein indicator**
- Dilute **nitric acid (HNO<sub>3</sub>)** to maintain ionic strength

#### Procedure:

1. **Pipette** a known volume of sodium chloride solution (e.g., 25 mL) into a conical flask.
2. Add **1–2 drops** of eosin Y or fluorescein indicator.
3. Add a few drops of **dilute HNO<sub>3</sub>** to maintain pH and ionic strength.
4. Titrate with **0.1 N AgNO<sub>3</sub>** from a burette while swirling the solution.
5. The solution will remain colorless during most of the titration.
6. At the **end point**, the color will change to **pale pink** or **red** due to adsorption of indicator onto excess Ag<sup>+</sup>-covered AgCl precipitate.

#### Reaction:



#### Calculation:

From the volume of AgNO<sub>3</sub> used, the amount of NaCl is calculated using:

$$1 \text{ mL of } 0.1 \text{ N AgNO}_3 \equiv 0.005844 \text{ g of NaCl}$$

### Advantages of Fajans Method:

- **Sharp and visible end point**
- No need for filtration or heating
- Suitable for **low concentrations of halides**

### Limitations:

- The **color change is sometimes subtle**
- **Photodecomposition** of AgCl may cause error
- Not suitable in **colored or turbid solutions**


## Complexometric Titration

### Definition of Complexometric Titration

Complexometric titration is a type of volumetric titration based on the formation of a stable **complex** between a **metal ion** and a **complexing agent (ligand)**. The most widely used reagent is **EDTA (ethylenediaminetetraacetic acid)**, a hexadentate ligand that forms stable chelates with most metal ions.

### Basis of Reaction

The titration involves the reaction:



Where **M<sup>n+</sup>** is a metal ion and **EDTA<sup>4-</sup>** is the complexing agent.

### Classification of Complexometric Titrations

Complexometric titrations can be classified into the following types based on the method of performing the titration:

#### 1. Direct Titration

- **Definition:** Metal ion is directly titrated with standard EDTA solution.
- **Conditions:** pH must be suitably buffered to ensure complete complex formation.
- **Example:** Estimation of calcium or magnesium in water using Eriochrome Black T indicator.

## 2. Back Titration

- **Definition:** A known excess of EDTA is added to the metal ion solution; the unreacted EDTA is titrated with a standard solution of another metal ion (e.g.,  $\text{Mg}^{2+}$  or  $\text{Zn}^{2+}$ ).
- **Use:** Applied when the direct titration is not feasible due to slow complex formation or precipitation of metal ions.
- **Example:** Estimation of aluminium or bismuth.

## 3. Replacement (Displacement) Titration

- **Definition:** A less stable metal-EDTA complex is displaced by a more stable metal ion.
- **Process:** EDTA is first added to form a weak complex with the analyte; then a stronger metal ion is added, which displaces the analyte from the complex. The displaced analyte is titrated.
- **Example:** Estimation of  $\text{Hg}^{2+}$  by replacing it with  $\text{Mg}^{2+}$  and titrating the released  $\text{Mg}^{2+}$ .

## 4. Indirect Titration

- **Definition:** The substance to be estimated is first converted into a form that can react with EDTA, and the metal ion content is determined indirectly.
- **Example:** Estimation of phosphate by precipitating it as  $\text{MgNH}_4\text{PO}_4$ , then determining the  $\text{Mg}^{2+}$  by EDTA titration.

## 5. Masking and Demasking Titrations

- **Masking:** A process of **preventing a particular ion** from reacting with EDTA by converting it into a non-reactive form.
- **Demasking:** The process of **releasing a masked ion** so it can be titrated.
- **Use:** Useful in selective titration of one metal ion in the presence of others.
- **Example:** Masking of aluminium ions using triethanolamine in the presence of calcium and magnesium.


# Metal ion indicators

### Definition:


Metal ion indicators are **organic dyes** that form **weakly dissociable, colored complexes** with metal ions. In complexometric titrations, they help **visually indicate the endpoint** by undergoing a distinct **color change** when the metal ion is completely complexed by the titrant (usually EDTA).

**Principle:**

At the start of the titration, the metal ion ( $Mn^+$ ) forms a weak, colored complex with the indicator (In):



During titration, EDTA (a stronger ligand) displaces the indicator from the metal ion due to its higher stability constant:



The **color change** of the indicator from metal-complex form to free indicator signals the **endpoint** of the titration.

**Ideal Properties of a Metal Ion Indicator:**

1. Forms a **weak and reversible complex** with metal ions.
2. Exhibits a **distinct color difference** between the free and complexed form.
3. Must not form a more stable complex than EDTA.
4. The pH range of the indicator must match the titration conditions.
5. Should be highly sensitive to even small changes in metal ion concentration.

**Common Metal Ion Indicators:**

Indicator	Metal Ion Detected	Color Change (Approx.)	pH Range
Eriochrome Black T	$Ca^{2+}$ , $Mg^{2+}$ (e.g., in water)	Wine red (complex) $\rightarrow$ Blue (free)	8–10
Murexide	$Ca^{2+}$	Pink (complex) $\rightarrow$ Purple (free)	11
Calcon	$Ca^{2+}$	Red (complex) $\rightarrow$ Blue (free)	12
Xylenol Orange	$Fe^{3+}$ , $Bi^{3+}$ , $Th^{4+}$	Red-orange (complex) $\rightarrow$ Yellow (free)	1–3
PAN (1-(2-Pyridylazo)-2-naphthol)	$Zn^{2+}$ , $Cu^{2+}$ , $Co^{2+}$	Red (complex) $\rightarrow$ Yellow (free)	3–5
Patton-Reeder Indicator	$Ca^{2+}$ in presence of $Mg^{2+}$	Red (complex) $\rightarrow$ Blue (free)	12



In complexometric titrations, especially when using EDTA to estimate metal ions, **interference** from other metal ions can lead to **inaccurate results**. To overcome this, specific reagents called **masking and demasking agents** are used.

## 1. Masking Agents

### Definition:

Masking agents are chemical substances that **selectively react with interfering metal ions** to form **stable, unreactive complexes**, thereby **preventing** them from reacting with the titrant (EDTA).

### Purpose:

To **"hide"** **interfering metal ions** without removing them physically from the solution, allowing accurate titration of the target metal ion.

### Mechanism:

The masking agent forms a complex with the interfering metal ion that is **more stable** than the one formed with EDTA, hence **inhibiting reaction** between the interfering metal ion and EDTA.

### Examples of Masking Agents:

Masking Agent	Masked Metal Ion(s)	Remarks
Potassium cyanide (KCN)	$\text{Cu}^{2+}$ , $\text{Cd}^{2+}$ , $\text{Zn}^{2+}$ , $\text{Ni}^{2+}$	Forms stable cyanide complexes
Triethanolamine (TEA)	$\text{Al}^{3+}$ , $\text{Fe}^{3+}$	Used in presence of $\text{Ca}^{2+}$ and $\text{Mg}^{2+}$
Sodium thiosulfate ( $\text{Na}_2\text{S}_2\text{O}_3$ )	$\text{Ag}^+$ , $\text{Hg}^{2+}$	Forms soluble thiosulfate complexes
Citric acid	$\text{Fe}^{3+}$ , $\text{Al}^{3+}$ , $\text{Mn}^{2+}$	Chelating ligand
Oxalate	$\text{Ca}^{2+}$ , $\text{Sr}^{2+}$	Forms insoluble calcium oxalate
Ascorbic acid	$\text{Fe}^{3+}$	Reduces to $\text{Fe}^{2+}$ which can be masked

## 2. Demasking Agents

### Definition:

Demasking agents are substances that **release a previously masked metal ion** from its masked complex, **restoring its ability** to react with EDTA or be titrated.

### Purpose:

To allow **sequential determination** of multiple metal ions in the same solution.

**Mechanism:**

The demasking reagent **breaks the complex** between the masking agent and metal ion, thereby freeing the metal ion for titration.

**Examples of Demasking Agents:**

Demasking Agent	Demasked Metal Ion(s)	Mechanism
Formaldehyde	$\text{Al}^{3+}$	Breaks Al–triethanolamine complex
Acetone	$\text{Fe}^{3+}$	Breaks Fe–citric acid complex
8-Hydroxyquinoline	$\text{Zn}^{2+}$ , $\text{Cu}^{2+}$	Precipitates metal ions
$\text{H}_2\text{S}$ (hydrogen sulfide)	$\text{Cd}^{2+}$ , $\text{Hg}^{2+}$	Precipitates as metal sulfides

**Applications of Masking and Demasking:**

1. Estimation of  $\text{Ca}^{2+}$  in presence of  $\text{Mg}^{2+}$  using **Patton-Reeder indicator** after masking  $\text{Mg}^{2+}$ .
2. Estimation of  $\text{Zn}^{2+}$  after masking interfering  $\text{Fe}^{3+}$  using **ascorbic acid**.
3. Sequential analysis of multiple metal ions in **pharmaceutical formulations** or **alloy samples**.

**1. Estimation of Magnesium Sulphate ( $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ )**

**Method Used:** Complexometric titration using **EDTA**

**Principle:** Magnesium ions ( $\text{Mg}^{2+}$ ) form a stable complex with EDTA in a buffered alkaline medium.


**Indicator:** Eriochrome Black T

**Buffer:** Ammonia-ammonium chloride buffer (pH 10)

---

**Procedure:**

1. Accurately weigh a quantity of magnesium sulphate and dissolve it in water.
2. Pipette an aliquot into a conical flask.
3. Add about 10 mL of pH 10 buffer solution.
4. Add 2–3 drops of **Eriochrome Black T** indicator; solution turns **wine red**.
5. Titrate with **0.01 M EDTA** until the color changes from wine red to **blue**.

**Reaction:**

**Calculation:**

1 mL of 0.01 M EDTA  $\equiv$  0.00247 g of  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$

**2. Estimation of Calcium Gluconate ( $\text{C}_{12}\text{H}_{22}\text{CaO}_{14} \cdot \text{H}_2\text{O}$ )**

**Method Used:** Complexometric titration using **EDTA**


**Principle:** Calcium ions ( $\text{Ca}^{2+}$ ) form a stable 1:1 complex with EDTA in alkaline medium.

**Indicator:** Murexide

**Buffer:** pH  $\sim$ 12 (NaOH used)

**Procedure:**

1. Dissolve the accurately weighed calcium gluconate in water.
2. Pipette an aliquot into a conical flask.
3. Add 1 mL of **NaOH** (to maintain pH  $\sim$ 12).
4. Add a small amount of **murexide indicator** (solution turns **pink**).
5. Titrate with **0.01 M EDTA** until the color changes from pink to **purple**.

**Reaction:****Calculation:**

1 mL of 0.01 M EDTA  $\equiv$  0.00427 g of calcium gluconate

**Precautions:**

- Maintain appropriate **pH** for each titration.
- Avoid contamination with other metal ions.
- Use freshly prepared indicator solutions.

**GRAVIMETRY****Definition:**

Gravimetry is a **quantitative analytical technique** in which the **amount of an analyte is determined based on the mass** of a solid product that is formed and weighed after chemical transformation and isolation.

**Principle of Gravimetric Analysis:**

The principle of gravimetry is based on the **measurement of mass**. It involves converting the analyte into a **pure, stable, and insoluble compound** of known composition, isolating it from the reaction mixture, drying (or igniting) it to a constant weight, and then **weighing** it accurately. The **amount of analyte** present in the original sample is calculated from the **mass of the precipitate** using **stoichiometry**.

### Key Conditions for Gravimetric Analysis:

1. The precipitate must be **pure, stable**, and of **known composition**.
2. The precipitation must be **complete** with minimum solubility of the product.
3. The analyte must be **completely and selectively precipitated**.
4. The product must be **easily filtered, washed, and dried/ignited**.
5. Precipitating agent must react **specifically** or **selectively**.

### Steps Involved :

#### 1. Sample Preparation

- The analyte (metal or ion to be estimated) is dissolved in a suitable solvent, usually water or dilute acid/base.
- If necessary, interferences are removed or masked before precipitation.

#### 2. Precipitation

- A **suitable precipitating agent** is added to the solution to form an **insoluble compound** of the analyte.
- Conditions such as **pH, temperature, and concentration** are adjusted to ensure complete and selective precipitation.
- The precipitate should form slowly to produce larger and purer crystals, which are easier to filter.

#### 3. Digestion

- The precipitate is allowed to stand in the hot mother liquor (solution from which it formed) for some time.
- This helps improve crystal size and purity by recrystallization (Ostwald ripening), making it easier to filter and wash.
- Digestion also helps in removing co-precipitated impurities.

#### 4. Filtration

- The precipitate is separated from the mother liquor using **filter paper** or a **porous crucible** (e.g., Gooch or sintered glass crucible).
- Care is taken to avoid loss of precipitate during transfer.

#### 5. Washing

- The precipitate is washed with cold distilled water or a suitable washing liquid to remove soluble impurities like electrolytes and mother liquor.
- The wash liquid should not dissolve or react with the precipitate.

## 6. Drying or Ignition

- The precipitate is either:
  - a) **Dried** at 100–120 °C to remove moisture (for hygroscopic or thermally unstable substances), or
  - b) **Ignited** at high temperatures (600–1000 °C) in a furnace to convert it into a stable form (e.g., metal oxide).
- This step ensures that the final product is in a known, constant composition suitable for weighing.

## 7. Weighing

- After drying or ignition, the crucible and precipitate are cooled in a **desiccator** to avoid moisture uptake.
- The final weight is accurately recorded using an **analytical balance**.

## 8. Calculation

- From the known stoichiometry of the precipitate and its mass, the **amount of analyte** in the original sample is calculated.

## Purity of the Precipitate: Co-precipitation and Post-precipitation

In gravimetric analysis, the **accuracy** of the result depends on obtaining a precipitate that is **pure and of known composition**. However, during precipitation, the product may contain **impurities** due to phenomena like **co-precipitation** and **post-precipitation**.

### 1. Co-precipitation

#### Definition:

Co-precipitation is the process by which **impurities are carried down along with the desired precipitate**, even though the impurity itself does not exceed its solubility limit under the experimental conditions.

#### Occurs:

- Simultaneously during the formation of the main precipitate
- Most significant in **initial stages of precipitation**

#### Causes of Co-precipitation:

- a) **Surface Adsorption** – Impurities get adsorbed on the surface of the precipitate crystals (e.g.,  $\text{Fe}^{3+}$  adsorbed on  $\text{BaSO}_4$ )
- b) **Occlusion** – Trapping of impurities inside growing crystals due to rapid precipitation
- c) **Inclusion** – Ions with similar size and charge as the primary ion occupy its position in the crystal lattice (e.g.,  $\text{K}^+$  replacing  $\text{NH}_4^+$ )

### Prevention/Minimization:

- Use **slow and controlled precipitation**
- Add precipitating reagent slowly with constant stirring
- Digest the precipitate in the mother liquor to promote crystal purity
- Wash the precipitate thoroughly

## 2. Post-precipitation

### Definition:

Post-precipitation occurs when a **second, unwanted precipitate forms after the desired precipitate** has already formed and started to settle. This new precipitate may contaminate the primary product.

### Occurs:

- **After** the main precipitation
- Due to **prolonged standing** of the mixture, allowing other ions to slowly form insoluble compounds

### Example:

- In estimation of calcium as calcium oxalate, **magnesium oxalate** may form if the solution is left standing too long

### Prevention/Minimization:

- **Filter the precipitate promptly** after digestion
- Avoid long standing time after precipitation
- Remove interfering ions beforehand (by masking or proper reagent selection)

### Key Differences Between Co-precipitation and Post-precipitation

Feature	Co-precipitation	Post-precipitation
Time of occurrence	During primary precipitation	After primary precipitation
Type of contamination	Impurities from same solution phase	A different compound precipitating later
Nature	Surface adsorption, inclusion, occlusion	Separate precipitation of another compound
Prevention	Slow precipitation, digestion, washing	Prompt filtration after digestion


## Estimation of Barium as Barium Sulphate (BaSO<sub>4</sub>)

(Gravimetric Method – As per B.Pharm PCI syllabus, Pharmaceutical Analysis)

---

### Principle:

Barium ions (Ba<sup>2+</sup>) are quantitatively precipitated from an aqueous solution using **dilute sulfuric acid (H<sub>2</sub>SO<sub>4</sub>)** as **insoluble barium sulfate (BaSO<sub>4</sub>)**. The BaSO<sub>4</sub> formed is a white crystalline precipitate which is filtered, washed, ignited to constant weight, and weighed. From the known stoichiometry, the amount of barium (or its salt) in the original sample is calculated.



### Requirements:

- Standard barium chloride solution (BaCl<sub>2</sub>)
  - Dilute H<sub>2</sub>SO<sub>4</sub>
  - Beakers, funnel, filter paper or sintered crucible
  - Desiccator and analytical balance
  - Hot water and wash bottle
- 

### Procedure:

1. Weigh accurately a sample containing barium (e.g., BaCl<sub>2</sub>·2H<sub>2</sub>O) and dissolve in distilled water in a beaker.
  2. Heat the solution to about 70–80°C.
  3. Add **dilute H<sub>2</sub>SO<sub>4</sub> slowly with constant stirring** to precipitate BaSO<sub>4</sub>.
  4. Digest the mixture (keep hot without boiling) for 30–60 minutes to promote crystal growth and reduce co-precipitation.
  5. Filter the precipitate through a **pre-weighed sintered glass crucible** or filter paper.
  6. Wash the precipitate with hot water to remove Cl<sup>–</sup> and SO<sub>4</sub><sup>2–</sup> impurities (test washings with AgNO<sub>3</sub> and BaCl<sub>2</sub>).
  7. Dry the crucible in an oven at 110–120°C or **ignite** in a muffle furnace at 600–800°C to constant weight.
  8. Cool in a **desiccator** and weigh.
-

### Calculations:

- Molar mass of  $\text{BaSO}_4$  = **233.39 g/mol**
- From the mass of  $\text{BaSO}_4$  obtained, calculate the amount of  $\text{Ba}^{2+}$  or  $\text{BaCl}_2$  present.

Amount of  $\text{Ba}^{2+} = (137.33/233.39) \times \text{Weight of BaSO}_4$

Amount of  $\text{BaCl}_2 \cdot 2\text{H}_2\text{O} = (244.26/233.39) \times \text{Weight of BaSO}_4$

### Precautions:


- Avoid excess  $\text{H}_2\text{SO}_4$  to prevent supersaturation and fine precipitate formation.
- Always wash until sulfate and chloride ions are absent in the filtrate.
- Ensure complete drying or ignition before weighing.
- Do not overheat the precipitate to prevent decomposition.

## Diazotisation Titration

### Basic Principles of Diazotisation Titration

Diazotisation titration is a **redox volumetric titration** method used primarily for the estimation of **primary aromatic amines**. In this method, the amine group ( $-\text{NH}_2$ ) is **converted into a diazonium salt** using **nitrous acid ( $\text{HNO}_2$ )** under cold acidic conditions.

The reaction proceeds as:



Where:

- $\text{Ar-NH}_2$  = Aromatic primary amine
- $\text{HNO}_2$  = Nitrous acid (generated in situ from sodium nitrite and HCl)
- $\text{Ar-N}_2^+ \text{Cl}^-$  = Diazonium chloride (stable at low temperature)

### Reagents Used

1. **Sodium nitrite ( $\text{NaNO}_2$ )** – Titrant
2. **Dilute hydrochloric acid ( $\text{HCl}$ )** – Acidic medium
3. **Potassium iodide and starch** – For iodometric end-point detection
4. **Primary aromatic amine solution** – Analyte

### Conditions for Diazotisation



- The reaction is performed at a **temperature of 0–5°C**, since diazonium salts are unstable at higher temperatures.
- Acidic medium (HCl) is essential for the formation of nitrous acid (HNO<sub>2</sub>) from sodium nitrite.



## Methods of Diazotisation Titration

### 1. Direct Method

- A known volume of the **aromatic primary amine** is taken in a flask.
- The solution is cooled in an ice bath (0–5°C).
- Sodium nitrite is titrated into the cold amine solution containing HCl.
- **End-point** is detected by adding starch-iodide paper or by external iodometric detection:
  - A drop of the reaction mixture is placed on starch-iodide paper.
  - The appearance of **blue color** (due to iodine-starch complex) indicates the end-point (excess nitrite oxidizes iodide to iodine).

### 2. Back Titration Method

- An excess of standard sodium nitrite solution is added to the cold acidified amine solution.
- After completion of reaction, the **excess nitrite** is determined by **iodometric titration** with **standard sodium thiosulphate**.



## Applications of Diazotisation Titration

- **Estimation of primary aromatic amines** like:
  - **Aniline**
  - **Sulphanilamide**
  - **Sulphadiazine**
  - **Sulphamethoxazole**
- Used in quality control of drugs containing aromatic amine groups.
- Official method in various pharmacopoeias for estimation of sulpha drugs.

## Advantages

- Highly **specific** for primary aromatic amines.

- Simple, sensitive, and does not require sophisticated instrumentation.
- Suitable for **bulk drug analysis** and routine quality control.

#### **Precautions**

- Strict temperature control (0–5°C) is essential to prevent decomposition of diazonium salts.
- Freshly prepared sodium nitrite solution should be used.
- Use accurate timing and indicators for end-point detection.